



A New Age of Healthcare



Since 1979

# Electro Surgical Unit Infusion Pump

HQ in Seoul  
Factory in Wonju

Transparent  
Ethical  
**Professional**  
Customer Satisfaction



---

Since 1979, Daiwha Corp. Ltd. has been a leader in the manufacturing and distribution of medical devices, Electro Surgical Unit and Infusion Pump.

In 2009 we enlarged our R&D center at our headquarters in Seoul and rebuilt our factory within DongWha Medical Device Complex in Wonju City. Through continuous internal R&D investments we have successfully stayed true to our core values while expanding into the competitive global market.

Using our company's philosophy for transparent, ethical, and professional management along with our social contribution mindset, our employees set the standard for excellence in customer satisfaction.

INFUSION SYSTEM  
**MEDIFUSION®**

04 DI-2X00 Series  
05 DI-4000 | DS-3000  
06 DI-6000 | DS-5000  
07 DI-6X00V Series

08 Accessories

SURGICAL SYSTEM  
**MEDITOM® | MEDISONIC®**

09 DT-150S/200S  
09 DT-300S/400S  
10 DE-400E/400  
11 DUB-100  
12 DU-Series

13 Accessories

DI-2000



DI-2200



DI-4000



DS-3000



DI-6000



DS-5000



DI-6000V



DI-6200V



CMS




STATION 1/2



IV SET



DROP SENSOR



DT-400S/300S



DT-200S/150S



DE-400E



DUB-100




DE-400



DU-143/137/127/115



DU-515/505



DU-343/337/327/315



DB-137



COMING  
UP NEXT





# VOLUMETRIC INFUSION PUMP DI-2000 | 2200 / MEDIFUSION®

DI-2000



DI-2200



DI-2000EN



DI-2000A



- Single & Dual Channel of Compact Design
- Double Door Lock System for Free Drop Protection
- IV Set Library | Feed Set Library (DI-2000EN Only)
- Event History Monitoring System
- Various Setting Modes
- Refurbished DI-2000A Animal Pump Available

VOLUMETRIC | SYRINGE INFUSION PUMP  
**DI-4000 | DS-3000** / MEDIFUSION®

DI-4000



DS-3000



Advanced Technology Infusion System  
Flexible Structure : Customized Station System  
Secondary Infusion Mode | Relay Infusion Mode  
Various Language Configuration  
Event History Monitoring System



VOLUMETRIC | SYRINGE INFUSION PUMP / MEDIFUSION®  
**DI-6000 | DS-5000**



DI-6000



DS-5000



- Convenience TFT-LCD & Jog Dial
- Docking Station System & Central Monitoring System
- Secondary Infusion Mode (DI-6000 Only) | Relay Infusion Mode
- Close Infusion System (Optional)
- Blood Transfusion (DI-6000 Only)
- PCA & TCI (DS-5000 Only)



VOLUMETRIC INFUSION PUMP  
**DI-6000V | 6200V** / MEDIFUSION®



DI-6000V



DI-6200V



- Convenience Touch Screen LCD
- Close Infusion System (Optional)
- Blood Transfusion
- Central Monitoring System (Wireless)
- PCA (DI-6000V Only)
- Automatically Controlled Brightness



## ▶ Central Monitoring System



▶ DAIWHA's NEW System facilitates convenient and safe for the medical staff and patient

- Real time monitoring of multiple infusion pumps at nurse station.

▶ Useful for checking multiple pump information

- Infusion/Syringe Pump information can be shown.
- The pumps are divided into four easy-to-see colors.
- The history is recorded according to the sequence of events and alarm.

▶ CMS with Products




## ▶ Station 1 & 2




- Both stand alone and up to 4 units docking station system are available.
- CMS at nurse station with docking station. (DI-6000 & DS-5000)

## ▶ Drop Sensor



- Detects free flow caused by abrupt device malfunctioning.
- Detects No Drop caused by an empty bottle or clogging of the tube.

## ▶ Infusion



- General IV Set (Standard)
  - ① Clip for DAIWHA Close Infusion System
- General IV Set (Needleless)
  - ② Needleless-type available
- Lightproof IV Set (Available)

ELECTRO SURGICAL UNIT  
**DT-400S | 300S / DT-200S | 150S** / **MEDITOM®**



DT-300S/400S



DT-200S/150S



- Practical & Convenient Interface
- Various Operation Modes
- Safety Function D.R.M System
- Automatic Patient Return Monitoring
- Highest Compatibility with Various Accessories

# ELECTRO SURGICAL UNIT DE-400E | DE-400 / MEDITOM®

PREMIUM

DE-400E



DE-400



▶ DAIWHA ADVANCED ENERGY DEVICES

MEDI-POLAR® DB-137



[Bipolar+Blade]

New Digital Technology  
Safety Function D.R.M System  
DAIWHHA APC Function & APS Mode  
Highest Compatibility with Various Accessories  
CUSA® Excel Handpieces Usable  
DAIWHHA's Advanced Energy System (DE-400 Only)

ULTRA ELECTROSURGICAL UNIT  
**DUB-100** / **MEDITOM®**



► DAIWHA ADVANCED ENERGY DEVICES

<p><b>MEDISONIC® DU-143/137/127/115</b></p>  <p>[Ultrasonic]</p>	<p><b>MEDISONIC® DU-343/337/327/315</b></p>  <p>[Ultrasonic]</p>	<p><b>MEDISONIC® DU-505/515</b></p>  <p>[Ultrasonic]</p>	<p><b>MEDIPOLAR® DB-137</b></p>  <p>[Bipolar+Blade]</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------

- New Digital Technology
- DAIWHA Complex Energy System
- DAIWHA APC Function & APS Mode
- Enhanced Seal Mode
- Usability of Daiwha Advanced Energy Device
  - MEDISONIC® Ultrasonic Energy Device
  - MEDIPOLAR® Bipolar Vessel Sealing Energy Device





SURGICAL SYSTEM  
**DU-Series** / MEDISONIC®

PREMIUM

**DU-1 Series**

[43 / 37 / 27 / 15 cm]



**DU-3 Series**

[43 / 37 / 27 / 15 cm]




**DU-5 Series**








[5 / 15 cm]



















Ultrasonic Energy Device  
User-Focused Design  
Enhanced Seal Mode  
Built-in Transducer  
Rapid Cutting & Lower Heating  
Different Shaft Lengths for User Needs

▶ ACCESSORIES





<p><b>1</b> </p> <p>· Disposable Twin Button Handle with Cable</p>	<p><b>2</b> </p> <p>· Reusable Handle with Cable</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------

<b>1</b>	<b>Ball Type</b> .....	
	· DEBS-07   EBS-07	
<b>2</b>	<b>Blade Type</b> .....	
	· DELS-06   DELS-07   ELS-07   DELS-10   DELS-17	
<b>3</b>	<b>Needle Type</b> .....	
	· DENS-07   ENS-07   DENS-10   DENS-17	
<b>4</b>	<b>Needle(Pin) Type</b> .....	
	· DEVNS-07   EVNS-07	
<b>5</b>	<b>Loop Type</b>	
	· DEWL-07   EWL-07 .....	
	· DEWL-13 (Semicircle_type) .....	
	· DEWL-13 (T_type) .....	

<p><b>1</b> </p> <p>· Disposable Dual Plate, Adult, REM, 3/5m Cable</p> <p>· Disposable Dual Plate, Infant, REM, 3/5m Cable</p> <p><b>3</b> </p> <p>· Disposable Dual Plate</p> <p><b>5</b> </p> <p>· Reusable Rubber Plate with Cable (CE/Audio JACK)</p>	<p><b>2</b> </p> <p>· Disposable Single Plate, Adult, NON-REM, 3/5m Cable</p> <p>· Disposable Single Plate, Infant, NON-REM, 3/5m Cable</p> <p><b>4</b> </p> <p>· Disposable Single Plate</p> <p><b>6</b> </p> <p>· Reusable Neutral Plate with Cable (CE/Audio JACK)</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p><b>1</b> </p> <p>· 310-100-05</p> <p><b>2</b> </p> <p>· 310-102-05</p> <p><b>3</b> </p> <p>· 310-110-05</p> <p><b>4</b> </p> <p>· 310-112-05   310-112-07</p> <p><b>5</b> </p> <p>· 310-115-07</p>	<p><b>6</b> </p> <p>· 310-120-10   310-130-10   310-130-20   310-140-10   310-150-10   310-260-10</p> <p><b>7</b> </p> <p>· 310-122-10   310-132-10   310-132-20   310-142-10   310-152-10</p> <p><b>8</b> </p> <p>· 310-160-10   310-160-20   310-170-05   310-170-10   310-180-05   310-180-10   310-190-05   310-190-10   310-200-10</p>	<p><b>9</b> </p> <p>· 310-162-10   310-162-20   310-172-10   310-182-10</p> <p><b>10</b> </p> <p>· 310-165-10   310-165-20   310-175-10   310-185-10   310-195-10   310-205-10</p>
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p><b>1</b> </p> <p>Double Footswitch (FS82)</p>	<p><b>2</b> </p> <p>Double Footswitch</p>	<p><b>3</b> </p> <p>Single Footswitch</p>
-------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------

<p><b>1</b> </p> <p>Reusable Bipolar Forceps Cable</p> <p><b>3</b> </p> <p>Snare Cable</p>	<p><b>2</b> </p> <p>Reusable Plate Cable</p> <p><b>4</b> </p> <p>Reusable Unipolar Handle with Cable</p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p><b>1</b> </p> <p>DT-200S, 150S DT-400S, 300S</p>	<p><b>2</b> </p> <p>DUB-100 DE-400</p>
----------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------



## CONTACT US

---

HOMEPAGE



YOUTUBE



KAKAO



INSTARGRAM

