



OSAKA DENTAL



Dental Unit

www.osakadental.com



OSA-A3

● ADVANTAGES

- 1) LED light in cabinet, more convenient for checking machine, other factories don't have / Panorama X-film viewer
- 2) 9 memories / 6 LED light / Large seat
- 3) Touch type instrument tray control system with LCD panel / Assistant control panel with 4-slot tool holder
- 4) Multifunction foot pedal / With Synchronized Movement

NEW



Touch screen instrument tray with self-test system (9 memories)
 Quartz watch and temperature display
 Three pressure meter for each hand-piece (patent design)
 24V voltage display & 10A DC current display
 With synchronized movement

触屏操作盘，带自检系统（9个记忆位）
 工具盘石英表&室温温度显示表
 三只手机气压独立调节显示表
 24V电压检测&10A直流监测显示
 带联动系统（靠背和座垫同时动，没有联动功能的只能动其中一个）



NEW



Lighting system in maintenance box
 箱体照明系统



Big Panorama X-film viewer
 全景观片灯



9bulbs color-adjusting (option) LED lamp 6bulbs color-adjusting
 9孔（标配）/6孔可调节色温双色LED灯



New multifunctional electric and air foot-pedal
 新型多功能电气一体脚踏



Large chair frame
 大椅架



2in 1 control box with USB connector
 治疗板程控器二合一
 USB接头总控，维修简单



Touch button assistant instrument tray Detachable cuspidor
 触摸感应按键辅控
 可手动拆卸的玻璃痰盂&注塑水龙头



Charging plug and USB connector for electrical machine
 器械盘220V充电插头&USB接口，方便小设备充电



International standard power cord in junction box
 品字形国际通用标准总电源线

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Operation lamp	6 LED
1 SET	Big panoramic X-film viewer	24V
DOUBLE	Armrest	

1 SET	Adjustable position	with 9 memories
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Foot pedal	multifunction
YES	Sanitary seamless cushion	
YES	Assistant tray with control button	Touch type
1 SET	Dentist stool	



OSA-4C

● ADVANTAGES

- 1) The most economic type
- 2) The chair frame is double rust thickened and stronger
- 3) Top sale in Middest East



Option:
Backrest control system for option
可选靠背



Removable
Ceramic spittoon option
可拆卸旋转陶瓷痰盂
可选玻璃痰盂



Self-water
supply button
自助供水系统



Halogen operation lamp
口腔冷光灯



PU cushion
Pu皮垫，便于清洁

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Halogen operation lamp	
1 SET	X-film viewer	24V
DOUBLE	Armrest	

1 SET	Adjustable headrest	
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Foot pedal	simple
YES	Sanitary seamless cushion	
NO	Assistant tray	
WEAK	Suction	No strong suction
1 SET	Dentist stool	



OSA-1-2018 (LED LAMP/ HALOGEN LAMP)

ADVANTAGES

- 1) Computer controlled 24VD, C.motor
- 2) The chair frame is double rust , thickened and stronger
- 3) X-film with sensor, other factories don't have
- 4) Top sale in Latin America

OPTIONAL/选配



Halogen operation lamp (complement)

口腔冷光灯/标配

HALOGEN LAMP



Sensor LED lamp

可调色温感应LED冷光灯, 手柄可拆卸, 便于消毒

LED LAMP



Removable and rotatable glass cuspidor
Ceramic spittoon option
可拆卸旋转玻璃痰盂
可选陶瓷痰盂



5-Seat Handpiece holder
五连挂架盒



Assistant operation tray
辅助控制工具盘



Self-water supply button
自助供水系统



Foot control behind
圆形控制脚踏



Big instrument tray
简洁方形大工具盘



Option: Backrest control system for option
可选靠背控制系统

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Halogen Operation lamp	
1 SET	X-film viewer	24V
DOUBLE	Armrest	

1 SET	Adjustable headrest	
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Round foot pedal	simple
YES	Sanitary seamless cushion	
YES	Assistant tray	
1 SET	Dentist stool	



OSA-A2

● ADVANTAGES

- 1) LED light in cabinet, more convenient for checking machine, other factories don't have
- 2) 9 memories / 6 LED light / Bigger seat
- 3) X-film with sensor, other factories don't have / Multifunction foot pedal
- 4) Assistant control panel with 4-slot tool holder / Touch type instrument tray control system



NEW



Lighting system in maintenance box
箱体照明系统



Detachable painting cuspidor is convenient for cleaning, Touch button assistant instrument tray
易于清洗/可拆卸可喷痰盂
触摸按键辅控



Separated electric and water system, orderliness tube convenient maintenance
水电分离, 箱体管排列
整齐易于维修

NEW



Touch button instrument tray(9 memories) pressure meter display
触摸按键器械盘9个记忆键
气压表



X-ray film viewer with sensor
手感应观片灯



New multifunction electric and air foot-pedal
新型多功能电气一体脚踏



6 bulbs adjustable LED lamp
6孔可调节色温LED灯



Bigger size chair frame
中等椅架

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Operation lamp	6 LED
1 SET	X-film viewer	24V, with sensor
DOUBLE	Armrest	

1 SET	Adjustable position	with 9 memories
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Foot pedal	multifunction
YES	Sanitary seamless cushion	
YES	Assistant tray with control button	Touch Type
1 SET	Dentist stool	



OSA-A1

● ADVANTAGES

- 1) LED light in cabinet, more convenient for checking machine, other factories don't have
- 2) 3 memories / Bigger seat / 6 LED light, if 4 LED light, reduce \$20 please
- 3) X-film with sensor, other factories don't have / Multifunction foot pedal
- 4) Assistant control panel with 4-slot tool holder



Lighting system in maintenance box
箱体照明系统



Assistant control box with control box system
时尚辅助盘控制系统



Multifunctional electric and air foot-pedal
多功能电气控制脚踏



Bigger size chair frame
中等椅架



New design maintain box
新型设计维修箱



Film view with sensor
豪华流线型铸铝平衡臂



6 bulbs adjustable LED lamp
6孔可调节色温LED灯



Big instrument tray with 3 memories
简洁方形工具盘 (带3个记忆位)

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Operation lamp	4 LED
1 SET	X-film viewer	24V, with sensor
DOUBLE	Armrest	

1 SET	Adjustable position	with 3 memories
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Foot pedal	multifunction
YES	Sanitary seamless cushion	
YES	Assistant tray with control button	
1 SET	Dentist stool	



OSA-1

● ADVANTAGES

- 1) The most top sale chair
- 2) Strong metal base with painting
- 3) Durable and stable quality



Cabinet with Painting Not only plastic injection, more durable
箱体是喷漆的，更耐用



Detachable comfortable armrest
可拆卸舒适款扶手



Sensor halogen lamp and key switch
带开关的卤素感应灯



Detachable injection cuspidor for convenient cleaning
可拆卸痰盂，注塑底座，方便清洁



Backrest control system is convenient for dentist.
便于医生和病人操作的靠背

SPECIFICATIONS

2 SETS	High-speed handpiece tubing with connector	2 holes/4 holes
1 SET	Low-speed handpiece tubing with connector	2 holes/4 holes
1 SET	3-way syringe	cold and hot water
1 SET	Water suction	
1 SET	Saliva ejector	
1 SET	Halogen operation lamp	
1 SET	X-film viewer	24V
SINGLE	Double armrest	

1 SET	Adjustable headrest	
1 SET	Water purified and heating system	
1 SET	Movable glass/ceramic cuspidor	
1 SET	Round foot pedal	simple
YES	Sanitary seamless cushion	
YES	Assistant tray	
1 SET	Dentist stool	



OSA-A6800



LED dental light
LED口腔灯



Physiological monitoring system
病人生理监护系统



Deluxe removable cleaning ceramic spittoon
Induction water supply system
Cuspidor impact protection system
可拆卸清洗陶瓷痰盂盘
感应供水系统
口水盆放碰撞保护装置



Digitized touch screen control system
Adjustable assistance arm and there are 3 positions
Dismountable holder and easy to clean
数字化触摸式按键控制系统
配合三段式可调高低转臂，操作更灵活
可拆卸器械支架，清洗灭菌更方便



Digitized touch screen control system
Digitized control system for ratio flow control of air/water
Hidden air lock for control balance arm system
数字化触摸式按键控制系统
数字化水气流量控制技术
可拆卸的挂架支架及气刹拉手，清洗灭菌更方便

Optional Color



蓝色 Blue



苹果绿 Apple Green



黄色 Yellow



酒红色 Claret

OSA-1-98A

全电脑连体式牙科综合治疗机
COMPUTER
CONTROLLED
INTEGRAL DENTAL UNIT

COMPUTER CONTROLLED INTEGRAL DENTAL UNIT



可旋转整体陶瓷痰盂 Rotatable Ceramic Spittoon



主控工具盘 Main Table

Technical Specification

- Supply voltage: AC110V/220V 60Hz/50Hz
- Motor voltage: DC24V
- Air pressure: 0.5MPa - 0.8MPa
- Water pressure: 0.2MPa - 0.4MPa
- Handpiece Tubing: 2/4-hole optional, Fiber optic tubing with extra charge

SPECIFICATIONS

- | | |
|-------|--|
| 1 SET | Thermostatic water heating system (applied to cup filler and three-way syringe) |
| 1 SET | Spittoon flushing and cup filling system (water supply amount can be preset) |
| 1 SET | Quiet electric lift chair (with adjustable pillow, single armrest) |
| 1 SET | Main table integrated with control panel, with 5-slot tool holder: 3 for handpiece tubings, 1 for three-way syringe, 1 open for optional use |
| 1 SET | LED film viewer |
| 1 SET | Assistant control panel with 5-slot tool holder: 2 for saliva suction/strong&weak, 1 for three-way syringe, 2 open for optional use |
| 1 SET | Rotatable ceramic spittoon |
| 1 SET | Brightness adjustable dental light (strong/weak switch control) |
| 1 SET | Multifunction pedal switch (for water/air supply, chair movement etc.) |
| 1 SET | Pure water bottle 600ml (removable, air pressurized, switch control) |
| 1 SET | Dentist stool |

Optional Equipment



牙科手机套装
Handpiece Set



内置洁牙机
Built-in Scaler



内置LED光固化机
Built-in LED Curing Light



口腔内窥镜一体机
Intra-oral Camera Set



医用空压机
Medical Compressor

NEW

OSA-1-28A -New

全电脑连体式牙科综合治疗机

COMPUTER
CONTROLLED
INTEGRAL DENTAL UNIT

COMPUTER CONTROLLED INTEGRAL DENTAL UNIT

多功能带冲痰、供水脚踏开关 " 只有您想不到的惊喜, 没有我们做不了的方案 ", 水气一体开关 : 方便、简捷、安全的设计理念。从平稳、高贵的连体落地式机箱到时尚、典雅、舒适的沙发真皮坐垫为您呈现科技与艺术的完美结合。



小推车款
Trolley Type



可旋转陶瓷痰盂及内置纸巾盒
Rotatable Spittoon & Built-in Tissue Box



高档触控屏 Deluxe Touch Screen



无影口腔灯 Shadowless LED Lamp (options)

SPECIFICATIONS - Side box unit is ground based type, high class inner tubing

- 1 SET Thermostatic water heating system (applied to cup filler and three-way syringe)
- 1 SET Spittoon flushing and cup filling system (supported by Italian solenoid valves, water supply amount can be preset)
- 1 SET Quiet electric lift folding chair (supported by High Class motors, large chair size, with adjustable pillow, double armrests)
- 1 SET Main table integrate with Touch screen control panel, with 5-slot tool holder: 3 for handpiece tubings, 1 for three-way syringe, 1 open for optional use
- 1 SET Chair movement system with 12 memory positions
- 1 SET Panoramic LED film viewer
- 1 SET Assistant control panel with 5-slot tool holder: 2 for saliva suction/strong&weak, 1 for three-way syringe, 2 open for optional use
- 1 SET Rotatable ceramic spittoon
- 1 SET Shadowless LED sensor dental light (cold light, gradual brightness adjustment)
- 1 SET Multifunction pedal switch (for water/air supply, chair movement etc.)
- 1 SET Double pure water bottles, 600ml each (removable, air pressurized, switch control)
- 1 SET Outer ground box for electric, air, water and drainage connection
- 1 SET Small glass tray and built-in tissue box
- 1 SET Deluxe dentist stool

Optional Equipment



牙科手机套装
Handpiece Set



内置洁牙机
Built-in Scaler



内置LED光固化机
Built-in LED Curing Light



口腔内窥镜一体机
Intra-oral Camera Set



医用空压机
Medical Compressor

OSA-1-68B

全电脑连体式牙科综合治疗机

COMPUTER
CONTROLLED
INTEGRAL DENTAL UNIT

COMPUTER CONTROLLED INTEGRAL DENTAL UNIT

现代化是它的代名词. 从可转动机箱、多功能脚踏、到水气电体总开关、手机水独立控制, 完美地体现了造型与功能的精致结合。



Technical Specification

- Supply voltage: AC110V/220V 60Hz/50Hz
- Motor voltage: DC24V
- Air pressure: 0.5MPa – 0.8MPa
- Water pressure: 0.2MPa – 0.4MPa
- Handpiece Tubing: 2/4-hole optional, Fiber optic tubing with extra charge



可旋转磨砂玻璃痰盂 Rotatable Mat Glass Spittoon



可旋转侧箱 Rotatable Side Box



特色工具盘 Specially Designed Table



豪华真皮皮垫 High Class Leather Upholstery

SPECIFICATIONS - Side box unit 45°rotatable.

- 1 SET Thermostatic water heating system (applied to cup filler and three-way syringe)
- 1 SET Spittoon flushing and cup filling system (water supply amount can be preset)
- 1 SET Quiet electric lift chair (with adjustable pillow, single armrest)
- 1 SET Main table integrated with control panel, with 5-slot tool holder: 3 for handpiece tubings, 1 for three-way syringe, 1 open for optional use
- 1 SET LED film viewer
- 1 SET Assistant control panel with 5-slot tool holder: 2 for saliva suction/strong&weak, 1 for three-way syringe, 2 open for optional use
- 1 SET Rotatable mat glass spittoon
- 1 SET Reflector LED sensor dental light (cold light, gradual brightness adjustment):
- 1 SET Multifunction pedal switch (for water/air supply, chair movement etc.)
- 1 SET Pure water bottle 600ml (removable, air pressurized, switch control)
- 1 SET Small glass tray
- 1 SET Dentist stool with backrest

Optional Equipment



牙科手机套装
Handpiece Set



内置洁牙机
Built-in Scaler



内置LED光固化机
Built-in LED Curing Light



口腔内窥镜一体机
Intra-oral Camera Set



医用空压机
Medical Compressor

color chart



Accessories for Dental Unit

Sealing machine	Sealing machine	whitening accelerator	Amalgamator	water Distiller	Monitor and oral camera	Sterilizer	Sputum aspirator	Steam sterilizer 18L
								
Microscope	Dental Ultrasonic Cleaner 5.0L	Handpiece	Dental Ultrasonic Cleaner 2.5L	Dental Ultrasonic Cleaner	Air Compressor for 6 (90-120L)	Air Compressor for 4 (60L)	Air Compressor for 1 (32L)	Suction Unit for 2
								