

OSAKADENT MEDICAL LIMITED


OSAKADENTAL

OSAKADENT MEDICAL LIMITED


Medium Green
 绿色、中等粒度



Item No.	DG009	DG010	DG010B	DG011	DG014	DG015	DG016	DG017	DG003A
Size (mm)	4 x 4	4 x 7	4 x 11	3 x 7	2.5 x 4.5	3.5 x 10	2.5 x 7	3.5 x 11	6.5 x 2



Item No.	DG001	DG002	DG002B	DG012	DG012B	DG013	DG004A	DG005A	DG004	DG005	DG006	DG007
Size (mm)	4 x 13	5 x 13	6.5 x 13	13 x 2	13 x 3.5	13 x 2	16 x 2/0.7	16 x 2	22 x 2/0.7	22 x 2	22 x 4.5	15 x 3.5


Extra Fine Yellow
 Gold shank

Suitable to adjusting Lithium Disilicate, All Ceramics, Zirconia of Occlusal, Lingual, Edge, Intracoronar, Proximal etc. Dentists will feel very convenient if they can prepare a few pieces beside the chairs.

专门针对诊所的全瓷、玻璃陶瓷打磨抛光、氧化锆牙调咬合、发育沟、窝沟、边缘、内冠、邻面、舌侧等。黄色、特细粒度。粒度细、不绷瓷。



Item No.	DY009EF	DY010EF	DY010BEF	DY011EF	DY014EF	DY015EF	DY016EF	DY017EF	DY003AEF	DY001EF	DY002EF	DY012EF	DY013EF
Size (mm)	4 x 4	4 x 7	4 x 11	3 x 7	2.5 x 4.5	3.5 x 10	2.5 x 7	3.5 x 11	6.5 x 2	4 x 13	5 x 13	13 x 2	13 x 2

OSAKADENT MEDICAL LIMITED

DuCoBur® 橡胶结合剂金刚石磨头
豆酷巴 (全瓷、全锆抛光磨头)

Diamond Polishers for Zirconia or AllCeramics
(Rubber Pregnant Diamond Polishing Burs)



Widely Used for Zirconia, Lithium Disilicate, Ceramics etc.
COMFORTABLE at the minimal pressure;
HEATLESS no need water cooling;
MORE DURABLE under our suggested PRM

不产热、无需水冷 更耐磨、经济实惠 手感好、很服帖 抛就亮、更轻松 型号齐、用途广、啥都抛得亮

Super Coarse
Black

黑色、特粗粒度



Item No.

D1S

D2AS

D2S

D3S

D3BS

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2

Coarse
Green

绿色、粗糙粒度



Item No.

D1C

D2AC

D2C

D3C

D3BC

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2

Medium
Blue

蓝色、中等粒度



Item No.

D1M

D2AM

D2M

D3M

D3BM

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2



D201M

5 x 16



D203M

6 x 16



D205M

12 x 2



D205BM

12 x 4



D207M

12 x 2.5




D13M

14 x 2



D17M

17 x 2.2



D26M

26 x 2.2

OSAKADENT MEDICAL LIMITED

 Fine Red
 红色、较细粒度



Item No.	D1F	D2AF	D2F	D3F	D3BF
Size (mm)	26 x 2	20 x 2	26 x 2	22 x 3	26 x 2



Item No.	D201F	D203F	D205F	D205BF	D207F	D13F	D17F	D26F
Size (mm)	5 x 16	6 x 16	12 x 2	12 x 4	12 x 2.5	14 x 2	17 x 2.2	26 x 2.2

 Extra Fine Yellow
 黄色、特细粒度



Item No.	D1EF	D2AEF	D2EF	D3EF	D3BEF
Size (mm)	26 x 2	20 x 2	26 x 2	22 x 3	26 x 2



Item No.	D201EF	D203EF	D205EF	D205BEF	D207EF	D13EF	D17EF	D26EF
Size (mm)	5 x 16	6 x 16	12 x 2	12 x 4	12 x 2.5	14 x 2	17 x 2.2	26 x 2.2

OSAKADENT MEDICAL LIMITED

 Ultra Fine
 Grey
 灰白色、特细粒度
 1



Item No.	D1UF	D2AUF	D2UF	D3UF	D3BUF
Size (mm)	26 x 2	20 x 2	26 x 2	22 x 3	26 x 2



D201UF	D203UF	D205UF	D205BUF	D207UF	D13UF	D17UF	D26UF
5 x 16	6 x 16	12 x 2	12 x 4	12 x 2,5	14 x 2	17 x 2,2	26 x 2,2

OSAKADENT