

ISO 9001 : 2000 CERTIFIED

ASCO

S I N C E : 1 9 4 8

Autoclaves

Sterilizers

Hospital Holloware

Hospital Utensils



Contents	Page No
Autoclave/Pressure Steam Sterilizer	
Autoclave/Pressure Steam Sterilizer, Aluminium - Pressure Cooker Type, Non Electric	1
Autoclave/Pressure Steam Sterilizer, Aluminium - Pressure Cooker Type, Electric	1
Autoclave/Pressure Steam Sterilizer Aluminium - Non Electric, Wingnut Type	2
Autoclave/Pressure Steam Sterilizer Aluminium - Electric, Wingnut Type	2
Autoclave/Pressure Steam Sterilizer, Stainless Steel	3
Autoclave/Pressure Steam Sterilizer, Vertical	3
Sterilizers	
Instrument Sterilizer, Stainless Steel, Non Electric	4
Instrument Sterilizer, Stainless Steel, Electric	4
Instrument Sterilizer, Stainless Steel, Electric, with Thermostat & Tray Lifting Arrangement	5
Bowl and utensil Sterilizer, Stainless Steel, Electric	5
Accessories - Autoclaves/Sterilizers	
Kerosene Stove, Single Burner	6
Kerosene Stove, Four Burners	6
Heating plate, single	7
Heating plate, four	7
Basket (for inserting inside Autoclave)	7
Immersion Heater	7
Syringe Rack (For Inserting Inside Autoclave)	7
Hospital Holloware	
Dressing Drums	8
Bowls Lotion, ss	9
Wash Basins, ss	9
Instrument tray with cover, ss	10
Catheter Tray with cover, ss	10
Tray Shallow, without cover ss	10
Kidney Tray, ss	10
Baby Tray, ss	11
Instrument Box, ss	11
Syringe and Instrument Box, Aluminium	11
Surgical Box, Aluminium	11
Instrument Box Deluxe, Aluminium	12
Midwifery Box, Aluminium	12
Bucket with Lid, ss	12
Water filter with candles, ss	12
Jars Forceps, ss	13
Jars Dressing/Ointment, ss	13
Jars- Dressing with Hinges, ss	13
Douche Can (Irrigator), ss	13
Gallipot with cover, ss	14
Sputum Mug with cover, ss	14
Feeding cup, ss	14
Measure, ss	14
Bed Pans & Urine Pots	
Bed Pan, Slipper Type - with cover, ss	15
Bed Pan, Perfection Type - with cover, ss	15
Bed Pan, Perfection Type seamless - with cover, ss	15
Bed Pan, Perfection Type seamless - with cover, ss	15
Fracture Bed Pan, Deluxe with cover, ss	16
Urinal/Urine Pot- Female, Deluxe , ss	16
Urinal/Urine Pot- Male, ss	17
Urinal/Urine Pot- Female, ss	17
Plastic Holloware	
Sputum Mug, Round, Plastic	18
Sputum Mug, Square, Plastic	18
Douche Set, Plastic	18
Kidney Tray, Plastic	18
Urinal/Urine Pot, Male, Plastic	19
Urinal/Urine Pot, Female, Plastic	19
Bed Pan, Slipper Type, Plastic	19
Bed Pan, Perfection Type, Plastic	19
Hospital Kitchenware	
Terms & Conditions inside back cover	20



Autoclave/Pressure Steam Sterilizer Aluminium - Pressure Cooker Type, Non Electric

Features:

- Aluminum deep drawn Seamless Construction
- Sturdy, Safe, Rust proof, Portable, Polished Shiny Finish
- For high pressure steam sterilization of medical materials such as dressings and surgical instruments.
- Max Pressure 15 psi/1.5 bar, Max temp 121° C
- Equipped with steam release/ control valve, high quality heat resistant Bakelite handles
- Lid Sealing by high pressure rubber gasket
- Pressure Gauge with color coding to display sterilization zones, graduated in kg/cm² and lb/inch²
- Equipped with Safety Pressure release weight valve
- Non-Electric

Sizes:

Code	Capacity (Litres) including Lid	Dia x Height (mm)	Dia x Height (inch)
AU2100	11	220 x 280	9 x 11
AU2200	15	300 x 225	12 x 9
AU2300	21	300 x 300	12 x 12
AU2400	24	300 x 350	12 x 14
AU2500	28	300 x 400	12 x 16
AU2600	38	350 x 400	14 x 16
AU2700	39	300 x 550	12 x 22

Autoclave/Pressure Steam Sterilizer Aluminium - Electric, Pressure Cooker Type

Features:

- Aluminum deep drawn Seamless Construction
- Sturdy, Safe, Rust proof, Portable
- For high pressure steam sterilization of medical materials such as dressings and surgical instruments.
- Max Pressure 15 psi/1.5 bar, Max temp 121° C
- Equipped with steam release/ control valve and high quality heat resistant Bakelite handles
- Lid Sealing by high pressure rubber gasket
- Pressure Gauge with color coding to display sterilization zones, graduated in kg/cm² and lb/inch²
- Equipped with Safety Pressure release weight valve
- Electric , 220v/50Hz,single phase

Size:

Code	Capacity (Litres) including Lid	Dia x Height (mm)	Dia x Height (inch)
AU3100	11	220 x 280	9 x 11
AU3200	15	300 x 225	12 x 9
AU3300	21	300 x 300	12 x 12
AU3400	24	300 x 350	12 x 14
AU3500	28	300 x 400	12 x 16
AU3600	38	350 x 400	14 x 16
AU3700	39	300 x 550	12 x 22



**Autoclave/Pressure Steam Sterilizer
Aluminium - Non Electric, Wignut/Knob Type**

Features:

- Aluminum deep drawn Seamless Construction
- Sturdy, Safe, Rust proof, Portable, Polished Shiny Finish
- Locking of lid by 6 nuts with high quality rubber gasket.
- For high pressure steam sterilization of medical materials such as dressings and surgical instruments.
- Max Pressure 15 psi/1.5 bar, Max temp 121° C
- Equipped with steam release/ control valve and high quality heat resistant Bakelite handles
- Pressure Gauge with color coding to display sterilization zones, graduated in kg/cm² and lb/inch²
- Equipped with Safety Pressure release weight valve
- Non-Electric

Sizes:

Code	Capacity (Litres) including Lid	Dia x Height (mm)	Dia x Height (inch)
AU4100	17	300 x 250	12 x 10
AU4200	21	300 x 300	12 x 12
AU4300	27	300 x 380	12 x 15
AU4400	39	350 x 380	14 x 15



**Autoclave/Pressure Steam Sterilizer
Aluminium - Electric, Wignut/Knob Type**

Features:

- Aluminum deep drawn Seamless Construction
- Sturdy, Safe, Rust proof, Portable
- Locking of lid by 6 nuts with high quality rubber gasket.
- For high pressure steam sterilization of medical materials such as dressings and surgical instruments.
- Max Pressure 15 psi/1.5 bar, Max temp 121° C
- Equipped with steam release/ control valve and high quality heat resistant Bakelite handles
- Pressure Gauge with color coding to display sterilization zones, graduated in kg/cm² and lb/inch²
- Equipped with Safety Pressure release weight valve
- Electric, 220v/50Hz, single phase

Sizes:

Code	Capacity (Litres) including Lid	Dia x Height (mm)	Dia x Height (inch)
AU5100	17	300 x 250	12 x 10
AU5200	21	300 x 300	12 x 12
AU5300	27	300 x 380	12 x 15
AU5400	39	350 x 380	14 x 15





Autoclave/Pressure Steam Sterilizer, Stainless Steel

Features:

- Stainless Steel deep drawn Seamless Construction
- Sturdy, Safe, Rust proof, Portable
- For high pressure steam sterilization of medical materials such as dressings and surgical instruments.
- Max Pressure 15 psi/1.5 bar, Max temp 121° C
- Equipped with steam release/ control valve and high quality heat resistant Bakelite handles
- Lid Sealing by high pressure rubber gasket
- Pressure Gauge with color coding to display sterilization zones, graduated in kg/cm² and lb/inch²
- Equipped with Safety Pressure release weight valve

Sizes:

Code	Capacity (Litres) including Lid	Dia x Height (mm)	Dia x Height (inch)	Power Source
AU6100	30	350 x 300	14 x 12	Electric 220v/50Hz
AU6200	30	350 x 300	14 x 12	Non Electric

Autoclave/Pressure Steam Sterilizer, Vertical

- Inner and outer chamber made of heavy gauge stainless steel
- Triple wall with steam jacket
- Jacket is insulated by high grade glass wool to minimise the temperature loss
- Fitted with safety valve, water level indicator, drain valve, pressure release valve
- Pressure adjustable from 15 to 20 psi
- Lid Sealing by high pressure jointless rubber gasket
- Electric 220v/50 Hz

Optional :

- Automatic Low water cut-off
- 0-60 min Mechanical timer
- Digital Temperature indicator
- Dial type analog temperature gauge
- Mechanical timer with alarm system

Model No.	Capacity in Litres	Chamber Size Dia x Height (mm)
AU451	35	300 x 500
AU452	52	350 x 550
AU453	75	400 x 600
AU454	175	550 x 750



Instrument Sterilizer, Non Electric, ss**Features:**

- Construction of Stainless steel
- For Steam Sterilization/disinfection of surgical materials /instruments
- With Handles for Lid and Base
- Available in various sizes, custom sizes can be made on request
- Works on external heat source such as Gas Burner/Kerosene Stove/Wood/Charcoal etc

Item Code	Size (mm)	Size (inch)
ST2100	200 x 125 x 75	8 x 5 x 3
ST2101	250 x 125 x 75	10 x 5 x 3
ST2102	300 x 150 x 100	12 x 6 x 4
ST2103	350 x 150 x 100	14 x 6 x 4
ST2104	400 x 150 x 100	16 x 6 x 4
ST2105	400 x 200 x 150	16 x 8 x 6
ST2106	450 x 200 x 150	18 x 8 x 6
ST2107	500 x 200 x 150	20 x 8 x 6
ST2108	500 x 250 x 175	20 x 10 x 7
ST2109	600 x 300 x 150	24 x 8 x 6
ST2110	600 x 300 x 150	24 x 12 x 6
ST2111	600 x 300 x 300	24 x 12 x 12

**Instrument Sterilizer, Electric, ss****Features:**

- Construction of Stainless steel
- For Steam Sterilization/disinfection of surgical materials /instruments
- With Handles for Lid and Base
- Available in various sizes, custom sizes can be made on request
- Electrically operated with immersion heater, Power source 220v/50hz
- Optional : Thermostat to control water temperature (please add T in the code)

Item Code	Size (mm)	Size (inch)
ST3100	200 x 125 x 75	8 x 5 x 3
ST3101	250 x 125 x 75	10 x 5 x 3
ST3102	300 x 150 x 100	12 x 6 x 4
ST3103	350 x 150 x 100	14 x 6 x 4
ST3104	400 x 150 x 100	16 x 6 x 4
ST3105	400 x 200 x 150	16 x 8 x 6
ST3106	450 x 200 x 150	18 x 8 x 6
ST3107	500 x 200 x 150	20 x 8 x 6
ST3108	500 x 250 x 175	20 x 10 x 7
ST3109	600 x 300 x 150	24 x 8 x 6
ST3110	600 x 300 x 150	24 x 12 x 6
ST3111	600 x 300 x 300	24 x 12 x 12

Seamless construction sterilizer

ST4100	300 x 150 x 125	12 x 6 x 5
ST4101	430 x 200 x 150	17 x 8 x 6
ST4102	500 x 200 x 150	20 x 8 x 6



Instrument Sterilizer, Stainless Steel, Electric, with Thermostat & Tray Lifting Arrangement

Features:

- Construction of Stainless steel
- For Steam Sterilization/disinfection of surgical materials /instruments
- With Handles for Lid and Base
- Available in various sizes, custom sizes can be made on request
- Electrically operated with immersion heater, Power source 220v/50hz

Item Code	Size (mm)	Size (inch)
ST5100	300 x 150 x 100	12 x 6 x 4
ST5101	350 x 150 x 100	14 x 6 x 4
ST5102	400 x 150 x 100	16 x 6 x 4
ST5103	400 x 200 x 150	16 x 8 x 6
ST5104	450 x 200 x 150	18 x 8 x 6
ST5105	500 x 200 x 150	20 x 8 x 6
ST5106	500 x 250 x 175	20 x 10 x 7
ST5107	600 x 300 x 150	24 x 8 x 6
ST5108	600 x 300 x 150	24 x 12 x 6
ST5109	600 x 300 x 300	24 x 12 x 12



Bowl and utensil Sterilizer, Stainless Steel, Electric

Features:

- Construction of Stainless steel
- For Steam Sterilization/disinfection of utensils, bowls etc
- Low water cut off device
- Fitted with Thermostat
- Available in various sizes, custom sizes can be made on request
- Electrically operated with immersion heater, Power source 220v/50hz

Item Code	Size (mm)	Size (inch)
ST6100	500 x 400 x 400	20 x 16 x 16
ST6101	600 x 400 x 400	24 x 16 x 16
ST6102	600 x 500 x 500	24 x 20 x 20
ST6103	600 x 600 x 600	24 x 24 x 24



Item Code : ASA101

Kerosene Stove - Single Burner

- Brass Tank
- Pressure type
- Tank capacity : 1.2 Litres approx
- Heavy Duty



Item Code : ASA102

Kerosene Stove - Four Burners

- Brass Tank
- Pressure type
- Tank capacity : 2.5 Litres approx
- Heavy Duty





Item Code : ASA103
Heating Plate, Single

- Electric Heating plate
- Fitted with coil type heating element
- Finished with attractive paint
- Power source: 200v/50hz



Item Code : ASA104
Heating Plate, Four

- Electric Heating plate
- Four plates in one unit
- Fitted with coil type heating element
- Finished with attractive paint
- Power source: 200v/50hz



Item Code : ASA105
Basket (for inserting inside Autoclave)

- Made of Aluminium
- Perforated for sterilization
- With Handle
- Available in various sizes

Item Code : ASA106
Immersion Heater

- For electric autoclaves and sterilizers



Item Code : ASA107
Syringe Rack (For Inserting Inside Autoclave)

- Made of Aluminium
- For placing of syringes inside autoclave
- Holds 36 syringes (4 syringes 5 ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles
- With Lid

Dressing Drums**Features:**

- Construction of Stainless steel
- Sterilizing Drums used for sterilizing and storing various types of dressings, instruments and surgical requisites
- Made of Stainless Steel
- Sets of perforation on drum body
- Slotted belt with clamp for loosening or tightening of belt to open or close the perforations on the body
- The clip when clamped will keep the belt in fair tension
- Available in various sizes, custom sizes can be made on request



Item Code	Dia × Height (mm)	Dia × Height (mm)
DD8100	150 x 100	6 x 4
DD8101	150 x 150	6 x 6
DD8102	200 x 150	8 x 6
DD8103	200 x 200	8 x 8
DD8104	225 x 112	9 x 4.5
DD8105	225 x 150	9 x 6
DD8106	225 x 225	9 x 9
DD8107	250 x 150	10 x 6
DD8108	250 x 250	10 x 10
DD8109	275 x 112	11 x 4.5
DD8110	275 x 225	11 x 9
DD8111	300 x 150	12 x 6
DD8112	300 x 250	12 x 10
DD8113	300 x 300	12 x 12
DD8114	325 x 150	13 x 6
DD8115	350 x 250	14 x 10
DD8116	375 x 225	15 x 9
DD8117	375 x 300	15 x 12
DD8118	375 x 150	15 x 6



Bowls Lotion, ss

Item Code	Dia in mm (approx)	Dia in inch (approx)	Capacity (ml)
HW1101	100	4	200
HW1102	125	5	380
HW1103	150	6	650
HW1104	175	7	1200
HW1105	200	8	1800

Wash Basins, ss

Item Code	Dia in mm (approx)	Dia in inch (approx)	Capacity (ml)
HW1201	300	12	3000
HW1202	325	13	4000
HW1203	350	14	5000
HW1204	375	15	6000
HW1205	400	16	7000



Instrument Tray with Cover, ss

Item Code	Size (Inch)	Size (mm)
HW2101	8x3x1.5	200x79x40
HW2102	8x6x2	200x150x50
HW2103	8.85x5x2	225x125x50
HW2104	9x6x2	225x150x50
HW2105	10x8x2	250x200x50
HW2106	11x7x2	275x175x50
HW2107	12x8x2	300x200x50
HW2108	12x10x2	300x250x50
HW2109	14x10x2	350x250x50
HW2110	15x12x2	375x300x50
HW2111	18x12x2	450x300x50
HW2112	8.85x3.28x1.65	222x82x41
HW2113	12.20x7.68x2.50	310x195x63

•Without cover available on request



Catheter Tray with cover, ss

Item Code	Size (Inch)	Size (mm)
HW2201	222x82x41	8½x3¼x1½
HW2202	430x100x70	17x4x2¾



Trays - Shallow, ss

Item Code	Size (Inch)	Size (mm)
HW2301	19x13x0.75	480x330x19
HW2302	13.50x9.75x0.625	343x247x16



Kidney Tray, ss

Item Code	Size (Inch)	Size (mm)	Capacity in ml (approx)
HW2401	6	150	225
HW2402	8	200	475
HW2403	10	250	825
HW2404	12	300	1000



Baby Tray, ss

Item Code	Dia in mm (approx)	Dia in inch (approx)
HW2501	450x30x75	18x12x3

Instrument Box, ss

Item Code	L×W×H (mm)
HW3101	180x80x30
HW3102	180x80x40
HW3103	180x80x50
HW3104	200x90x30
HW3104	200x90x40
HW3105	200x90x50
HW3106	250x100x40
HW3107	250x100x50
HW3108	250x100x60



Syringe and Instrument Box, Aluminum

Item Code	L×W×H (mm)
HW3201	80x35x20
HW3202	100x50x25
HW3202	120x55x25

Surgical Box, Aluminum

Item Code	L×W×H (mm)
HW3301	305x216x150
HW3302	465x285x100



Jars - Forceps, ss

Item Code	Dia×H×Base (mm)	Dia×H×Base (inch)
HW4301	50×125×100	2×5×4
HW4302	50×175×100	2×7×4

**Jars - Dressing/Ointment, ss**

Item Code	Dia × H (mm) approx	Capacity (Ltr.)
HW4401	100×130	1.00
HW4402	110×150	1.50
HW4403	165×125	2.00
HW4404	135×170	2.50
HW4405	160×200	4.00

**Jars - Dressing with Hinges, ss**

Item Code	Dia × H (mm) approx	Capacity (Ltr.)
HW4406	150×150	2.50
HW4407	200×150	4.50

**Douche Can (Irrigator), ss**

Item Code	Capacity (Ltr.) approx
HW4501	1.50
HW4502	2.00

Optional : Douche tube & Douche fitting set can be supplied on request



Bucket with Lid, ss

Item Code	Capacity (Ltr.) approx
HW4101	10
HW4102	12
HW4103	15
HW4104	18
HW4104	22

• without lid also available on request

Instrument Box Perforated Deluxe, Aluminum

Item Code	L × W × H (mm)	L × W × H (inch)
HW3401	350 × 200 × 80	14 × 8 × 3
HW3402	485 × 285 × 100	19.4 × 11.4 × 4



Water Filter with Candles, ss

Item Code	Capacity in liters (both containers combined)	No. of filter candles
HW4201	15	1
HW4202	18	1
HW4203	21	2
HW4204	24	2
HW4204	27	3
HW4205	29	3

Midwifery Box, Aluminum

Item Code	L × W × H (mm)
HW3501	305 × 216 × 150



Gallipot with Cover, ss

Item Code	Capacity (oz)
HW4601	4
HW4602	6
HW4603	8



Sputum Mug with Cover, ss

Item Code	Dia × H (mm)	Dia × H (inch)
HW4701	100 x 50	4 x 2



**Item Code : HW 4801
Feeding Cup, ss**



Measure, ss

Item Code	Capacity (ml)
HW4901	500

BEDPANS & URINE POTS

ISO 9001 : 2000 CERTIFIED

asco

SINCE : 1948



Item Code : HW1501

Bed Pan, Slipper Type - with cover, ss

Item Code : HW5102

**Bed Pan, Perfection Type -
with cover, ss**



Item Code : HW5103

Bed Pan, Perfection Type, Seamless - with cover, ss



Item Code : HW5104

**Bed Pan, Round Type -
with cover, ss**



Item Code : HW5104
**Fracture Bed Pan, Deluxe -
with cover, ss**



Item Code : HW5201
Urinal / Urine Pot - Female, Deluxe, ss



Item Code : HW5202
Urinal / Urine Pot - Male, ss



Item Code : HW5203
Urinal / Urine Pot - Female, ss



Item Code : HW5204
Urinal / Urine Pot - Male, ss



Item Code : HW5205
Urinal / Urine Pot - Female, ss

Item Code : HW6201
Sputum Mug Round, Plastic



Item Code : HW6202
Sputum Mug Square, Plastic

Item Code : HW6301
Douche Set, Plastic



Kidney Tray, Plastic

Item Code	Size (inch)	Size (mm)
HW6401	8	200
HW6402	10	250



Item Code : HW6501
Urinal/Urine Pot, Male, Plastic

Item Code : HW6502
Urinal/Urine Pot, Female, Plastic



Item Code : HW6601
Bed Pan, Slipper Type, Plastic

Item Code : HW6602
Bed Pan, Perfection Type, Plastic

Item Code : HW6603
Bed Pan, Perfection Type, Plastic



HW6602



HW6603

Features

- High quality stainless steel
- Long lasting
- Elegant looks
- Easy to clean
- Well finished
- Available in various designs and sizes



- Serving Utensils and Bowls
- Ladles, Serving spoons
- Cooking Pans and Bowls
- Food trays, Serving Trays
- Spoon, forks, knives
- Storage Drums
- Jars, Mugs
- Tea pots, Kettles

“Apothecaries Sundries Manufacturing Company” (ASCO) an ISO 9001:2000 certified company was founded in year 1948 and has since been developing and marketing various Medical and Surgical equipment. With its zest to achieve latest innovations in technology and fulfilling wide variety of customer's requirement it has developed a wide range of Autoclaves, Sterilizers & Hospital Holloware to meet broad spectrum of customer's requirements.

The Company has already secured for itself a position of reliability and prominence as makers of high quality products of international Standards.

The Company has the required infrastructure and all the testing facilities to ensure that the best quality products are provided to the customer. While most of the products are CE certified, a few are also FDA (USA) Listed. Also the company holds WHO GMP Compliance, in addition to ISO 9001 certification. The company has been recognized as One Star Export House, based on its export performance by Government of India.

Worldwide patronage : There is practically no country of the world where the 'ASCO' products have not found way at one time or the other, either direct or through re-exports by our customers. At present, the 'ASCO' products are being exported to over forty countries.

Meeting individual requirements : Every order, howsoever big or small is complied within accordance with individual requirement.

Twin Passion : Commitment to quality as well as holding the priceline is the twin passion with 'ASCO'

Assets : Ever growing clientele and excellent reference-both bank & trade-are our assets we feel justly proud of.

TERMS & CONDITIONS

- All sizes, dimensions, values etc. in this catalogue are approximate. We reserve the right to upgrade/improve/alter the specifications of our products without notice.
- We reserve the right to supply different colour/design than those shown in the catalogue, which depends on the availability and trend of raw material at the time of supply.
- Use of all our products is restricted to professionals who are qualified and trained to use them. We are not responsible for any damages which may occur due to improper or unauthorized use. In case of malfunctioning of the product, the company's liability is limited to the maximum extent of the cost of equipment.
- It is users responsibility to keep the products and packing away from the reach of unauthorised people, children and animals and dispose them off in an environment friendly manner
- All transactions shall be governed by jurisdiction under Delhi courts.
- You must keep proper record of your sales. This will assist product re-calls (if required)
- All orders are accepted subject to existence of conditions that enable smooth compliance. No liability shall lie against us or any of our representatives, if any orders, remain uncomplied with or partially complied with or delayed compliance for reasons beyond our control.

International Quality Products



**WHO
GMP**
COMPLIANT

ISO
13485:2003
CERTIFIED

ISO
9001:2000
CERTIFIED

**FDA
USA**
LISTINGS



Government of India recognised

"STAR EXPORT HOUSE"

Apothecaries' Sundries Mfg. Pvt. Ltd

ASCO House, 229, Udyog Vihar Phase IV
Gurugram-122015, Haryana, India
www.ascomedical.com