

# CATALOGUE



— SWEDISH DENTAL SUPPLIES AB —



# INDEX

<b>DISPOSABLES</b>	<b>2</b>
<b>ENDODONTICS</b>	<b>5</b>
<b>CROWNS</b>	<b>6</b>
<b>GENERAL DENTISTRY</b>	<b>8</b>



**Postal address**

Södervägen 30  
SE-232 52 Åkarp, Sweden

**Phone**

+46 40 46 02 01

**E-mail**

[info@swedent.se](mailto:info@swedent.se)

**Website**

[www.swedent.se](http://www.swedent.se)



## ARTICULATING PAPER, BLUE/RED 65 MY

Double Check, with blue and red colour on opposite side of the paper and with finger grip.

In 3 shapes:

- No. 1010 - Horseshoe, Box with 6 books of 12 sheets
- No. 1011 - Straight, Box with 12 books of 12 sheets
- No. 1012 - Curved, Box with 12 books of 12 sheets



## PLASTIC TOOTH CLEANER

- \* Replace the traditional wooden tooth pick.
- \* Can be bent without breaking.
- \* Flexible, follows the tooth surface.
- \* Lightly teased surface helps clean between the teeth.
- \* Hygienical packed in dispenser of about 100 pcs.
- \* 10 dispensers in each display cart.


No. 1035 - 10 boxes x 100 pcs



## SALIVA EJECTORS SWEFLEX

Environmental friendly material - for single use only!  
 Hygienic suction body with superior suction capacity reducing aerosol in treatment room.  
 Anatomical design - keeps tongue in place. Very steady!  
 Preshaped and flexible - one size - fits all patients.  
 Conical connection  $\varnothing$  6-7mm. (No adaptors needed).  
 Environmental friendly material, Polyethylen, products of combustion only CO<sub>2</sub> and H<sub>2</sub>O. Hypoallergic.

No. 1342 - Bag of 100 pcs



## SALIVA EJECTORS BIO ORIGINAL HYGOFORMIC


With tongue guard in one sling unit. Adjustable in many ways. Each bag incl. 2 adaptors.  
 Made of green PE, which is a fossil-free bio based material made from suger canes.

Packed in cart. of 20 bags x 100 pcs.

No. 1315 - Bag of 100 pcs



# DISPOSABLES



## ASPIRATOR TUBES

Without angle. Length 130 mm.  
Packed in carton of 10 bags x 100 pcs.

No. 1330 - Bag of 100 pcs white




## BARMAN WOODEN WEDGES

Anatomically shaped and polished for a better and safer filling procedur. Easy to insert between the teeth. Made of splinter resistant birchwood.

No. 1384 - Box of 200 pcs, small

No. 1385 - Box of 200 pcs, large

No. 1386 - Box of 200 pcs, small & large



## ADAPTOR FOR SALIVA EJECTORS

For connection of saliva ejectors into ejection system.  
Ø 11m.

No. 1320 - Bag of 10 pcs




## ROOT CANAL PINS

A flexible, residual plastic burn-out pin with good retention in the impression materials. Colour: Pink and green.

No. 3290 - Box of 60 pcs pink


No. 3280 - Box of 60 pcs green



## HYGOSURGE

Surgical aspirator to fit directly to 11 mm system.  
Apex: 3,6 mm. Length: 165 mm. Color: White. Autoclavable

No. 1337 - Bag of 25 pcs




## MICROTIP

Autoclavable suction device for Hygosurge. A tiny blue autoclavable tip that turns the Hygosurge into a unique suction device for micro oral surgery and fine endodontic work.

Colour: Light-blue.


No 1338 - Bag of 25 pcs



## ARTICULATOR - OCCLUDATOR

Disposable of plastic. For inlays, crowns and small bridges. Also for fitting study models. Colour: Blue.

No. 1021 - Box of 25 pcs



## APPLICATOR BRUSHES

Disposable dental applicators for precision application of: Cavity lining - Etching - Sealants - Bonding - Adhesives - Varnish application - Haemostatic liquids - Disclosing agents

No. 1395 - Handle white, bag of 5 pcs

No. 1396 - Mini-brushes, bag of 100 pcs




## ANGLE WIDER AND LIP HOLDER

Easy to place and suitable for photography, oral surgery, orthodontic bonding and restorative work. Made of polycarbonate and sterilizable in autoclave at 134°C.

No. 4370 - Angle wider, adult. 2 pcs

No. 4371 - Lip holder, adult. 2 pcs



## DISPOSABLE PLASTIC SYRINGE

Save messy, time-consuming clean up. For impression materials, cements etc. Volume 5 cc.

No. 1387 - Box of 100 pcs



## DAPPEN-DISHES


Disposable dappen dishes withstand any liquid using in a dental operator, as fluorides, acrylics, chloroform or prophy paste. Ideal as prophy cups for mixing acrylics and amalgams. 4 different colours to avoid confusion and contamination of materials. The wide base remains steady in use. The side tap may be bent up to form handle for use as prophy dish.

No. 1390 - Bag of 50 pcs yellow

No. 1391 - Bag of 50 pcs blue

No. 1392 - Bag of 50 pcs white

No. 1393 - Bag of 50 pcs red



## PIPETTE

Stands any liquid in a dental office, makes it possible to effectively rinse the root canals. Long, thin, but very firm tip. Volume 0,5 ml.

No. 1042 - Box of 100 pcs, non sterilized





# ENDODONTICS



## RUBBER DAM

In 3 colours, **dark, light and green**. Withstand for autoclaving upto 118°C (224°F). Vacuum packed for prolonged shelf life.

### IN ROLLS:


- No. 3081 - Thin light, roll of 450 x 15 cm
- No. 3060 - Thin green, roll of 450 x 15 cm
- No. 3083 - Medium light, roll of 450 x 15 cm
- No. 3061 - Medium green, roll of 450 x 15 cm

## RUBBER DAM

In 3 colours, **dark, light and green**. Withstand for autoclaving upto 118°C (224°F). Vacuum packed for prolonged shelf life.

### IN CUT PIECES:

- No. 3086 - Thin light, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3070 - Thin green, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3088 - Medium light, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3071 - Medium green, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3091 - Thin light, 6" x 6", 15,0 x 15,0 cm, 36 pcs
- No. 3075 - Thin green, 6" x 6", 15,0 x 15,0 cm, 36 pcs
- No. 3093 - Medium light, 6" x 6", 15,0 x 15,0 cm, 36 pcs
- No. 3076 - Medium green, 6" x 6", 15,0 x 15,0 cm, 36 pcs



## RUBBER DAM FRAME N/O

Od plastic and radiolucent. No need to remove when X-raying. Autoclavable upto 118°C (224°F)

No. 3100



## RUBBER DAM CLAMPS

Numbers:	15	24	27
	18	25	28
	22	26	29
	30	138	201
	31	139	202
	51	200	203
	204	207	212
	205	208	0
	206	209	00

No. 3110 - S.S. White pattern. Available for incisors, bicuspid and molars



## RUBBER DAM CLAMPS HOLDER PALMER

No. 3112

## RUBBER DAM PUNCH AINSWORTH

No. 3111



# CROWNS



## TEMPORARY POLYCARBONATE CROWNS

The crowns are made of non-toxic polycarbonate. This plastic is due to its molecule structure extremely unsusceptible, even to the strain of important temperature changes and other physical or chemical influences. These very flexible and thin crowns can be cut, formed, bent and pressed without breaking or ripping.

They form stable compounds with all acrylates and ZnO-tooth cements. Further these polycarbonate crowns are superior because of their good cervical fit, precise anatomy and shade stability.

**No. 2190** - Anteriors and Premolars assortment box of 180 crowns, without mould guides, 3 crowns x 60 sizes

**No. 2195** - Anteriors and Premolars assortment box of 180 crowns, with mould guides, 3 crowns x 60 sizes

**No. 2198** - Anteriors and Premolars assortment box of 360 pcs crowns, with mould guides, 6 crowns x 60 sizes

**No. 2215** - Refill box of 5 crowns of any mould

**No. 2220** - Loose crowns in bulk, min. 100 pcs/each mould

## MOULD GUIDES

Upper Centrals		Upper Laterals		Cuspids		Lower Anteriors		Bicuspids	
Right	Left	Right	Left	Right	Left	Long	Short	Large	Small
100	103	2	200	300	303	60	69	40	50
101	102	20	29	301	302	61	68	41	51
10	19	21	28	30	39	62	67	42	52
11	18	22	27	31	38	63	66	43	53
12	17	23	26	32	37	64	65	44	54
13	16	24	25	33	36				
14	15			34	35				





## CROWN FORMS, strip crowns SYSTEM C

The strip crowns are made out of extremely thin, super elastic, form stable non-toxic material. Their excellent anatomical design and occlusion surface reduces the necessary correcting work-done with a scalpel or scissor - to a minimum. The long tooth-necks make optimal fitting possible. The clear material allow control of the hardening and the colour. They can also be used with UV-light hardening materials. The crown material does not bind with any of the used plastics, and thus the crowns are easy to remove. The crowns are available in assortment boxes of 60 respectively 120 pieces and in refill boxes. The crowns are available in all sizes of anteriors and premolars.

- No. 2140 - Anteriors and Premolars assortment box of 60 crowns, 1 crown x 60 sizes
- No. 2141 - Anteriors and Premolars assortment box of 120 crowns, 2 crowns x 60 sizes

- No. 2150 - Refill box of 3 crowns of any mould
- No. 2155 - Loose crowns in bulk, min. 100 pcs each mould

1 WIDTH: 9,75 mm  A-1 L CUS. XL	2 WIDTH: 8 mm  A-2 L LAT. XL	3 WIDTH: 10 mm  A-3 L CEN. XL	4 WIDTH: 10 mm  A-4 R CEN. XL	5 WIDTH: 8 mm  A-5 R LAT. XL	6 WIDTH: 9,75 mm  A-6 R CUS. XL	7 WIDTH: 5 mm  F-1 L CEN. UD	8 WIDTH: 10 mm  G-1 L CEN. XL	9 WIDTH: 10 mm  G-2 R CEN. XL	10 WIDTH: 10 mm  L-1 L CEN. XL	11 WIDTH: 10 mm  L-2 R CEN. XL	12 WIDTH: 7,1 mm  P-1 L-R BICUS.
13 WIDTH: 9 mm  B-1 L CUS. L	14 WIDTH: 7,5 mm  B-2 L LAT. L	15 WIDTH: 9,5 mm  B-3 L CEN. L	16 WIDTH: 9,5 mm  B-4 R CEN. L	17 WIDTH: 7,5 mm  B-5 R LAT. L	18 WIDTH: 9 mm  B-6 R CUS. L	19 WIDTH: 5 mm  F-2 R CEN. UD	20 WIDTH: 10 mm  H-1 L CEN. L	21 WIDTH: 10 mm  H-2 R CEN. L	22 WIDTH: 9,5 mm  M-1 L CEN. XL	23 WIDTH: 9,5 mm  M-2 R CEN. XL	24 WIDTH: 7 mm  P-2 L-R BICUS.
25 WIDTH: 8 mm  C-1 L CUS. M	26 WIDTH: 6,75 mm  C-2 L LAT. M	27 WIDTH: 8,75 mm  C-3 L CEN. M	28 WIDTH: 8,75 mm  C-4 R CEN. M	29 WIDTH: 6,75 mm  C-5 R LAT. M	30 WIDTH: 8 mm  C-6 R CUS. M	31 WIDTH: 4,5 mm  F-3 L-R CEN LD	32 WIDTH: 9 mm  I-1 L CEN. M	33 WIDTH: 9 mm  I-2 R CEN. M	34 WIDTH: 9,5 mm  N-1 L CUS. XL	35 WIDTH: 9,5 mm  N-2 R CUS. XL	36 WIDTH: 6,9 mm  P-3 L-R BICUS.
37 WIDTH: 7,5 mm  D-1 L CUS. MS	38 WIDTH: 6,25 mm  D-2 L LAT. MS	39 WIDTH: 8 mm  D-3 L CEN. MS	40 WIDTH: 8 mm  D-4 R CEN. MS	41 WIDTH: 6,25 mm  D-5 R LAT. MS	42 WIDTH: 7,5 mm  D-6 R CUS. MS	43 WIDTH: 6,75 mm  F-4 L-R BICUS. M	44 WIDTH: 8,5 mm  J-1 L CEN. M	45 WIDTH: 8,5 mm  J-2 R CEN. M	46 WIDTH: 9 mm  O-1 L CUS. XL	47 WIDTH: 9 mm  O-2 R CUS. XL	48 WIDTH: 6,8 mm  P-4 L-R BICUS.
49 WIDTH: 6,75 mm  E-1 L CUS. S	50 WIDTH: 5,75 mm  E-2 L LAT. S	51 WIDTH: 7,5 mm  E-3 L CEN. S	52 WIDTH: 7,5 mm  E-4 R CEN. S	53 WIDTH: 5,75 mm  E-5 R LAT. S	54 WIDTH: 6,75 mm  E-6 R CUS. S	55 WIDTH: 7,25 mm  F-5 L-R BICUS. L	56 WIDTH: 8 mm  K-1 L CEN. S	57 WIDTH: 8 mm  K-2 R CEN. S	58 WIDTH: 6,3 mm  P-7 L-R BICUS.	59 WIDTH: 6,5 mm  P-6 L-R BICUS.	60 WIDTH: 6,7 mm  P-5 L-R BICUS.



## SCREW POSTS GOLD PLATED

Conical shaped to avoid unwanted wedge effect. Electrolytically plated with a thick layer of 24K gold. In 20 sizes and 4 lengths.

Sizes:		
Short,	8,0 mm	No. 1-6
Medium,	9,5 mm	No. 1-6
Long	12,0 mm	No. 1-6
Extra Long,	14,0 mm	No. 4 and 6

- No. 4399 - Ass. box of 60 pcs (3 pcs of 20 sizes+key+wrench)
- No. 4400 - Ass. box of 120 pcs (6 pcs of 20 sizes+key+wrench)
- No. 4401 - Ass. box of 240 pcs (12 pcs of 20 sizes+key+wrench)
- No. 4402 - Refill box of 12 pcs

- No. 4420 - Cross key, separate
- No. 4421 - Hollow key, separate




## SCREW POSTS PURE TITANIUM

A tissue compatible. Corrosion resistant and no chemical reactions in the human body. In 20 sizes and 4 lengths.

Sizes:		
Short,	8,0 mm	No. 1-6
Medium,	9,5 mm	No. 1-6
Long,	12,0 mm	No. 1-6
Extra Long	14,0 mm	No. 4 and 6

- No. 4410 - Ass. box of 60 pcs (3 pcs of 20 sizes+key+wrench)
- No. 4411 - Ass. box of 120 pcs (6 pcs of 20 sizes+key+wrench)
- No. 4413 - Ass. box of 240 pcs (12 pcs of 20 sizes+key+wrench)
- No. 4412 - Refill box of 6 pcs



## TITANIUM DENTINE RETENTION PINS

Titanium is tissue compatible and has no corrosion or galvanic reaction associated with other metals.

- \* non corrosive titanium
- \* smooth-cutting
- \* no stress
- \* maximum retention
- \* soft-shearing
- \* cut from solid rod
- \* self-centering
- \* superior stability

The retention pins are supplied in handy 12-packs with one drill, in diameters 0,75 mm (medium) and 0,60 mm (small)

- No. 4450 - Kit box of 12 pcs + 1 drill diameters 0.75 mm
- No. 4451 - Kit box of 12 pcs + 1 drill diameters 0.60 mm



### CODE RINGS ROULETT 225 rings

- 9 colours x 25 rings. For colour coding hand instruments.
- Silicon Rubber rings, fit all hand instruments.
  - Sterilizable up to 240°C / 464°F
  - In 9 different colours for easy identification.

No. 4050 - Assortment box. 225 rings




### CODE RINGS

- For colour coding hand instruments.
- Silicon Rubber rings, fit all hand instruments
  - Sterilizable up to 240°C / 464°F
  - In 11 different colours for easy identification.

Box of 50 rings - one colour.

- |                   |                   |                   |
|-------------------|-------------------|-------------------|
| No. 4058 - Pink   | No. 4062 - Brown  | No. 4066 - White  |
| No. 4059 - Purple | No. 4063 - Green  | No. 4067 - Yellow |
| No. 4060 - Black  | No. 4064 - Grey   | No. 4068 - Red    |
| No. 4061 - Blue   | No. 4065 - Orange |                   |

No. 4070 - Assortment box of 80 rings = 10 rings each of 8 colours.



### TAPEN TELL

A surgical and dental colour coding tape in 8 colours. Spool of 3200 x 7 mm. Autoclavable.

- |                  |                   |
|------------------|-------------------|
| No. 4350 - Black | No. 4354 - Orange |
| No. 4351 - Blue  | No. 4355 - Red    |
| No. 4352 - Brown | No. 4356 - White  |
| No. 4353 - Green | No. 4357 - Yellow |




# GENERAL DENTISTRY



## PROTECTIVE GOGGLES, KLEERSITE

Made of Polycarbonate with side shield. Transparent colour.  
Packed in white box.

No. 4271



## COMPOSITE INSTRUMENT

Double plane plastic for all composite materials. In 2 sizes:  
Autoclavable upto 118°C (224°F)

No. 6120 - No. 1, White. No. 6121 - No. 2, Blue

Instruments especially suited for light curing composite.  
Autoclavable upto 118°C (224°F)

No. 6122 - No. 1, Clear. No. 6123 - No. 2, Clear



## NAPKIN HOLDER JUMBO

A coloured napkin holder with large grips which makes it easy to mount it on the napkin. Cold-sterilizable. Light-weight.

No. 4710 - Red  
No. 4711 - White  
No. 4713 - Blue  
No. 4714 - Yellow  
No. 4715 - Green

## NAPKIN HOLDER

No. 4700 - with white plastic string




## DAPPEN DISHES

In 5 colours. Packed in box of 36 pcs, one colour.

No. 4243 - Glass clear  
No. 4245 - Glass pink  
No. 4246 - Glass yellow  
No. 4247 - Glass green  
No. 4241 - Glass dark-blue






## FILM MOUNT CARDS PLASTIC TRANSLUCENT

**Model 6** - Cards for 6 films bite-wing display  
No. 4390 - 105 x 120 mm. Package of 100 pcs

**Model 14** - Card for 14 films for general mouth survey  
No. 4391 - 135 x 200 mm. Package of 100 pcs

**Model 16** - Card for 16 films for assembly of full mouth survey  
No. 4392 - 135 x 200 mm. Package of 100 pcs

**Model 16 BW** - Card for 16 films bite-wing display  
No. 4394 - 135 x 200 mm. Package of 100 pcs



## FILM HANGER

No. 4559 - For 1 film


No. 4560 - For 10 films

No. 4561 - For 14 films



## MATRIX RETAINER

No. 4285 - Toffelmire model. Universal



## IWANSON MEASURING DEVICE type

Made of stainless steel. Very useful for checking the wall thickness of gold and porcelain crowns.

Doublesided with an accuracy of 1/10 of a mm.

No. 4447






### COTTON PELLET DISPENSER


No. 4720 - Chromium plated brass

No. 4730 - Plastic dome, separate



### ARTICULATING PAPER FORCEPS

No. 4750 - Stainless steel




### STAINLESS INSTRUMENT TRAY

No. 5672 - 284 x 183 mm



### BUR STAND ROTO OF PLASTIC

No. 5525 - For 18 FG and 54 Latch burs with dom.



### BUR STAND

40 x 40 x 12 mm

No. 5445 - For 12 FG + 12 Latch burs




### PELLET CONTAINER

50 x 27 x 10 mm

No. 5710 - Stainless steel




# GENERAL DENTISTRY



### BUR BOX

40 x 40 mm

No. 5511 - Stainless steel



### BUR BOX


50 x 27 x 10 mm

No. 5510 - Stainless steel (standard)



### ENDO-CLEAN


No. 5447 - For cleaning files in use



### ENDO INSTRUMENT STAND

100 x 50 x 25 mm

No. 5444 - For 38 files



### ENDO INSTRUMENT STAND, MINI

50 x 50 x 25 mm


No. 5443 - For 20 files



### BUR BOX

40 x 75 mm


No. 5512 - Stainless steel



### BOX

With 3 partitions. 50 x 27 x 10 mm

No. 5501 - Stainless steel



### BOX FOR POINTS

With 2 partitions. 50 x 27 x 10 mm.

No. 5500 - Stainless steel



### MEDICAMENT CUP

50 x 24 x 10 mm.

No. 5690 - Stainless steel

SWEDISH DENTAL SUPPLIES AB IS SERVING DENTAL DISTRIBUTORS WORLD WIDE WITH HIGH QUALITY DENTAL PRODUCTS SINCE 1966



# — SWEDISH DENTAL SUPPLIES AB —

*Postal address*

Södervägen 30  
SE-232 52 Åkarp, Sweden

*Phone*

+46 40 46 02 01

*E-mail*

info@swedent.se

*Website*

www.swedent.se