

ASTM RATED DISPOSABLE GLOVES



Item Code	Size	Dimension (mm)	Cuff (mm)	Palm Thickness (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Packing
-----------	------	----------------	-----------	---------------------	---------------------	-----------------------	---------



550-CP	XS	78mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	100pcs/bx 10bxs/ctn
650-CP	S	84mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
750-CP	M	93mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
850-CP	L	101mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
950-CP	XL	110mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	

Color:

Meets ASTM D3578

Feature: Advance "STRETCH" formulation providing maximum comfort for extended period usage.

Utilization: Fit-feel disposable latex gloves for normal examination procedures.



550-CPF	XS	78mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	100pcs/bx 10bxs/ctn
650-CPF	S	84mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
750-CPF	M	93mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
850-CPF	L	101mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
950-CPF	XL	110mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	

Color:

Meets ASTM D3578

Feature: The only "TACKY GRIP" disposable latex gloves in market for best gripping and handling.

Utilization: Fit-feel disposable latex gloves for normal examination procedures.



550-C	XS	78mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	100pcs/bx 10bxs/ctn
650-C	S	84mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
750-C	M	93mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
850-C	L	101mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	
950-C	XL	110mm ± 3	min. 240	min. 0.13	min. 0.10	min. 0.15	

Color: (B), (Y), (G), (K), (A)

Meets ASTM D3578

Fragrance Available: Bubblegum, Grape, Lemon, Mint, Strawberry **FRAGRANCE™ is required**

Utilization: "TACKY GRIP" disposable latex gloves for normal examination procedures.



550-CPNF	XS	< 80mm ± 10	min. 240	min. 0.07	min. 0.06	min. 0.08	100pcs/bx 10bxs/ctn
650-CPNF	S	80mm ± 10	min. 240	min. 0.07	min. 0.06	min. 0.08	
750-CPNF	M	95mm ± 10	min. 240	min. 0.07	min. 0.06	min. 0.08	
850-CPNF	L	110mm ± 10	min. 240	min. 0.07	min. 0.06	min. 0.08	
950-CPNF	XL	≥ 110mm ± 10	min. 240	min. 0.07	min. 0.06	min. 0.08	

Color:

Meets ASTM D6319

Feature: Finger textured for perfect grip.

Utilization: Fit-feel disposable nitrile gloves for normal examination procedures.

ASTM RATED DISPOSABLE GLOVES



Item Code	Size	Dimension (mm)	Length (mm)	Palm Thickness (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Packing
-----------	------	----------------	-------------	---------------------	---------------------	-----------------------	---------



Latex Examination Glove - Powdered

550-ECOCP	XS	78mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	100pcs/bx 10bxs/ctn
650-ECOCP	S	84mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	
750-ECOCP	M	93mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	
850-ECOCP	L	101mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	

Color:

Meets ASTM D3578

Utilization: Fit-feel disposable latex gloves for normal examination procedures.



Latex Examination Glove - Powder Free

550-ECOCPF	XS	78mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	100pcs/bx 10bxs/ctn
650-ECOCPF	S	84mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	
750-ECOCPF	M	93mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	
850-ECOCPF	L	101mm ±3	min. 240	min. 0.11	min. 0.08	min. 0.12	

Color:

Meets ASTM D3578

Utilization: Fit-feel disposable latex gloves for normal examination procedures.



Double Duty Glove - Powder Free

0-311-55X	5.5	72mm ±4	min. 280				2pairs/pouch 25pouches/bx 8bxs/ctn
0-311-60X	6.0	77mm ±5	min. 280				
0-311-65X	6.5	83mm ±5	min. 280				
0-311-70X	7.0	89mm ±5	min. 280	Inner: 0.12 ±0.02	Inner: 0.12 ±0.02	Inner: 0.13 ±0.02	
0-311-75X	7.5	95mm ±5	min. 280	Outer: 0.13 ±0.02	Outer: 0.12 ±0.02	Outer: 0.15 ±0.02	
0-311-80X	8.0	102mm ±6	min. 280				
0-311-85X	8.5	108mm ±6	min. 280				
0-311-90X	9.0	114mm ±6	min. 280				

Color: (Inner), (Outer)

Meets ASTM D3577

Feature: Two layers of gloves per hand for "DOUBLE PROTECTION".

STERILE



High Risk Glove - Powder Free

550-HR	XS	76mm ±3	300mm ±10	0.20mm		0.22mm	50pcs/bx 10bxs/ctn
650-HR	S	84mm ±3	300mm ±10	0.20mm		0.22mm	
750-HR	M	94mm ±3	300mm ±10	0.20mm	N/A	0.22mm	
850-HR	L	105mm ±3	300mm ±10	0.20mm		0.22mm	
950-HR	XL	113mm ±3	300mm ±10	0.20mm		0.22mm	

Color:

Meets ASTM D3578

Utilization: Safety disposable latex gloves, extended length and extra thickness for high risk environment.

ASTM RATED DISPOSABLE GLOVES



Item Code	Size	Dimension (mm)	Length (mm)	Palm Thickness (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Packing
0-111-55X	5.5	72mm ± 4	min. 270	min. 0.18	min. 0.17	min. 0.20	50pairs/bx 6bxs/ctn
0-111-60X	6.0	77mm ± 5	min. 275	min. 0.18	min. 0.17	min. 0.20	
0-111-65X	6.5	83mm ± 5	min. 275	min. 0.18	min. 0.17	min. 0.20	
0-111-70X	7.0	89mm ± 5	min. 280	min. 0.18	min. 0.17	min. 0.20	
0-111-75X	7.5	95mm ± 5	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-111-80X	8.0	102mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-111-85X	8.5	108mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-111-90X	9.0	114mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	



Color: **STERILE R**
 Meets ASTM D3577
 Feature: Advance "STRETCH" formulation providing maximum comfort for extended period usage.
 Utilization: Disposable latex gloves for surgery procedures.



0-211-55X	5.5	72mm ± 4	min. 270	min. 0.18	min. 0.17	min. 0.20	50pairs/bx 6bxs/ctn
0-211-60X	6.0	77mm ± 5	min. 275	min. 0.18	min. 0.17	min. 0.20	
0-211-65X	6.5	83mm ± 5	min. 275	min. 0.18	min. 0.17	min. 0.20	
0-211-70X	7.0	89mm ± 5	min. 280	min. 0.18	min. 0.17	min. 0.20	
0-211-75X	7.5	95mm ± 5	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-211-80X	8.0	102mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-211-85X	8.5	108mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	
0-211-90X	9.0	114mm ± 6	min. 285	min. 0.18	min. 0.17	min. 0.20	

Color: **STERILE R**
 Meets ASTM D3577
 Feature: Advance "STRETCH" formulation providing maximum comfort for extended period usage.
 Utilization: Disposable latex gloves for surgery procedures.

NON ASTM RATED DISPOSABLE GLOVES



ASTM RATED DISPOSABLE FACE MASK



Item Code

Type

Inspirational (ISO)

Packing



ASTM LEVEL 1
FACE MASK

Pro-Q 200 3-ply Surgical Face Mask

2-400 Fitting: Ear Loop/Tie On
Material: Standard/Sensitive
User: Adult/Child

Low Fluid Resistance	80 mmHg
Filtration Efficiency	BFE ≥ 95% PFE ≥ 95% @ 0.1 micron
Breathability - Delta P	< 4.0 mm H ₂ O/cm ²
Flame Spread	Class 1

STANDARD:
50pcs/bx
20bxs/ctn
SENSITIVE:
40pcs/bx
20bxs/ctn

Color: ●●●●●●●●

Meets EN14683 Rating - Type I Standard

Usage: Ideal for procedures where low amounts of fluid, spray or aerosols are produced.

Recommendation: Choose "SENSITIVE" material for maximum skin comfort.



ASTM LEVEL 2
FACE MASK

Pro-Q 200 Plus 3-ply Surgical Face Mask

2-406 Fitting: Ear Loop/Tie On
Material: Standard/Sensitive
User: Adult/Child

Moderate Fluid Resistance	120 mmHg
Filtration Efficiency	BFE ≥ 98% PFE ≥ 98% @ 0.1 micron
Breathability - Delta P	< 5.0 mm H ₂ O/cm ²
Flame Spread	Class 1

STANDARD:
50pcs/bx
20bxs/ctn
SENSITIVE:
40pcs/bx
20bxs/ctn

Color: ●●●●●●●●

Meets EN14683 Rating - Type IIR Standard

Usage: Ideal for procedures where moderate to light amounts of fluid, spray or aerosols are produced.

Recommendation: Choose "SENSITIVE" material for maximum skin comfort.



ASTM LEVEL 3
FACE MASK

Pro-Q 200 Ultra 3-ply Surgical Face Mask

2-407 Fitting: Ear Loop/Tie On
Material: Standard/Sensitive
User: Adult/Child

High Fluid Resistance	160 mmHg
Filtration Efficiency	BFE ≥ 98% PFE ≥ 98% @ 0.1 micron
Breathability - Delta P	< 5.0 mm H ₂ O/cm ²
Flame Spread	Class 1

STANDARD:
50pcs/bx
20bxs/ctn
SENSITIVE:
40pcs/bx
20bxs/ctn

Color: ●●●●●●●●

Meets EN14683 Rating - Type IIR Standard

Usage: Ideal for procedures where heavy to moderate amounts of fluid, spray or aerosols are produced.

Recommendation: Choose "SENSITIVE" material for maximum skin comfort.



ASTM LEVEL 1
FACE MASK

3 Piece 3-ply Face Mask

2-433 Fitting: Ear Loop/Tie On
Material: Standard
User: Adult/Child

Low Fluid Resistance	80 mmHg
Filtration Efficiency	BFE ≥ 95% PFE ≥ 95%
Breathability - Delta P	< 4.0 mm H ₂ O/cm ²
Flame Spread	Class 1

3pcs/pouch
50pouches/bx

Color: ●●●●●●●●

Meets EN14683 Rating - Type I Standard

Usage: Ideal for procedures where low amounts of fluid, spray or aerosols are produced.

NON ASTM RATED DISPOSABLE FACE MASK



1-ply/2-ply Paper Face Mask



2-ply Non Woven Face Mask



Cone Face Mask



General Face Mask



Utility Face Mask

NON ASTM RATED FACE SHIELD



Protective Eyewear



Protective Goggles



Protective Shield



Disposable Face Shield

PROTECTION APPAREL

Category: Bodywear



SMS Surgical Gown



Reinforced SMS Surgical Gown



Isolation Gown



PE Coated Isolation Gown



Coverall Suit



PE Coated Coverall Suit



Plastic Apron LDPE



Plastic Apron HDPE



Laboratory Coat



Plastic Gown



Spallace Jacket

Category: Footwear



Standard Shoe Cover



Antiskid Shoe Cover



Boot Cover



Plastic Shoe Cover

PROTECTION APPAREL

Category: Headwear

Available in weight 14G/30G



Color: Blue Green Yellow Pink White

J-703/2-707

Round Cap



Color: Blue Green Yellow Pink White

J-702

Clip Cap

Available in Type ELASTIC/TIE



Color: Blue Green

J-705

Surgeon Cap



Color: Blue

J-704

Astronaut Cap



Color: White

J-706

Beard Cover



Color: Blue

J-500

Plastic Cap



Color: White

J-03/707

Hairnet



Color: White

J-06007

Shower Cap

Category: Handwear



Color: Blue

J-708

Sleeve Cover



Color: Blue

J-500

Plastic Sleeve Cover

STERILIZATION PRODUCTS



Size Available:
 57x133mm
 60x130mm
 70x160mm
 90x207mm
 90x260mm
 130x280mm
 130x283mm
 130x258mm
 190x260mm
 200x350mm
 300x474mm

Made of MEDICAL GRADE PAPER 70G
 suitable for STERILIZATION



00113/00130/00160/00207/00260/00280/00283/00258/00260/00350/00474

Sterilization Pouches

Size Available:
 40x40cm
 40x45cm
 50x50cm
 60x60cm
 70x70cm
 80x80cm
 100x100cm
 120x120cm
 120x140cm
 120x160cm
 120x180cm
 140x140cm

Made of MEDICAL GRADE PAPER 40G
 suitable for STERILIZATION



00400/00405/00500/00600/00700/00800/01000/01200/01200/01400/01200/01600/01200/01800/01400

SMS Wrap

Size Available:
 50mmx200m
 75mmx200m
 100mmx200m
 125mmx200m
 150mmx200m
 200mmx200m
 250mmx200m
 300mmx200m
 350mmx200m
 400mmx200m

Made of MEDICAL GRADE PAPER 70G
 suitable for STERILIZATION




00500/00750/01000/01250/01500/02000/02500/03000/03500/04000

Sterilization Flat Reel

Size Available:
 40x40cm
 60x60cm
 80x80cm
 70x70cm
 90x90cm
 100x100cm
 100x140cm
 120x120cm
 120x140cm
 120x160cm
 120x180cm

Made of 100% VIRGIN WOOD PULP MEDICAL GRADE PAPER 60G



Color:
 ● Blue
 ● Green

00400/00600/00800/00700/00900/01000/01000/01200/01200/01400/01200/01400/01200/01600/01200/01800

Medical Crepe Paper

Size Available:
 75mmx100m
 100mmx100m
 150mmx100m
 200mmx100m
 250mmx100m
 300mmx100m
 350mmx100m
 400mmx100m

Made of MEDICAL GRADE PAPER 70G
 suitable for STERILIZATION



00750/01000/01500/02000/02500/03000/03500/04000

Sterilization Gasketed Reel

Size Available:
 12mmx50m
 13mmx50m
 25mmx50m

Perfect indicator, easy removed and replaced with no adhesive residue



01250/01350/02550

Plasma Indicator Tape

Size Available:
 50mmx50m
 75mmx50m
 100mmx50m
 150mmx50m
 200mmx50m
 250mmx50m
 300mmx50m
 350mmx50m
 400mmx50m

Made of 100% MEDICAL GRADE PAPER
 suitable for PLASMA STERILIZATION



00500/00750/01000/01500/02000/02500/03000/03500/04000

Sterilization Plasma Reel

WOUND CARE PRODUCTS

Category: Bandage & Accessories

Available in
size 1"/2"/3"/4"/6"



3-105/3-107/3-109/3-110/3-110

W.O.W Bandage 4.5m

Available in
size 9cm/7.5cm/10cm/15cm



3-200/3-201/3-202/3-203

Crepe Bandage 4.5m

Available in
size 1"/2"/3"/4"/6"



3-204/3-205/3-206/3-207/3-208

Self-Seal Bandage

Available in
size 9cm/7.5cm/10cm/15cm/20cm



3-233/3-234/3-235/3-236/3-237

Orthopaedic Bandage 1.75m

Available in
size 9cm/7.5cm/10cm/15cm/20cm



3-209/3-209/3-209/3-209/3-209

POP Bandage 2.7m

Available in
size 90x90x127cm



3-210

Triangular Bandage

Available in
size 12.5cm/25cm/50cm



3-211/3-211/3-211

Surgical Tape

Available in
size 2.5cm/5cm/7.5cm/10cm/15cm



3-212/3-212/3-212/3-212/3-212

Zinc Oxide Adhesive Plaster 5m

Category: Gauze & Sponges

Available in
type WITH/WITHOUT
X-RAY DETECTABLE



3-403/3-403

BPC Gauze Roll 4-ply

Available in
size 5x5cm/7.5x7.5cm/10x10cm



Also available in 5 ply
and 10 ply
3-405/3-405/3-405

Gauze Sponge 8-ply

Available in
size 5x5cm/7.5x7.5cm/10x10cm



Also available in 5 ply
and 10 ply
3-406/3-406/3-406

Sterile Gauze Sponge 8-ply

Available in
size 1x1cm/2x2cm



3-407/3-407

Flat Gauze



3-408

Non Woven Sponge Roll 100m

Available in
size 5x5cm/7.5x7.5cm/10x10cm



3-409/3-409/3-409

Non Woven Sponge 4-ply

Available in
size 5x5cm/7.5x7.5cm/10x10cm



3-411/3-411/3-411

Sterile Non Woven Sponge 4-ply

Available in
size 10x10cm



3-412

Tetra Gauze Sponge

WOUND CARE PRODUCTS

Category: Cotton



Cotton Wool Roll 500g



Cotton Wool Roll 400g



Combined Dressing Roll 300g



Sterile Cotton Ball 0.5g



Sterile Cotton Ball 5g



Cotton Ball 0.5g



Cotton Swab



Cotton Tip Applicator

Category: Other Accessories



Surgical Tape Dispenser



Bandage Scissor

HOSPITAL PRODUCTS

Category: Hospital Disposable



Basic Dressing Set



Alcohol Swab



Identification Tag



Disposable Brief



Disposable Underpad



Disposable Pillow Case



Disposable Bedsheet



Disposable Blanket



Compressed Towel



Compressed Facial Towel

Category: Hospital Non Disposable



Tourniquets



Pen Light



Pill Box

DENTAL SUPPLIES & ACCESSORIES

Category: Dental Accessories



Dental Tweezer

S-700



Dental Probe

S-707



Dental Mirror

S-702



Dental Tray

S-700



Impression Tray

SSD-709



Cheek Retractor

SSD-707



Bite Tray

SSD-702



Fluoride Tray

SSD-709



Denture Box

A777



Retainer Box

A777



Mixing Bowl

A701



Disposable Bob Clip

A700-D



Dental Applicator

A700-A



Film Mount

A701



Mixing Pad

A701/A702





Dappen Dish

S-200

DENTAL SUPPLIES & ACCESSORIES



Category: Suction & Aspirator Tip



DENTAL SUPPLIES & ACCESSORIES

Category: Endo & Impression Material



Intra Oral Tip



Mixing Tip



Mixing Gun

Category: Barrier Sleeve & Cover



Barrier Film



Air/Water Syringes Tip Cover



Air/Water Syringes Sleeve



Handpiece Sleeve



T-Style Lamp Sleeve




L-Style Lamp Sleeve



X-Ray Sleeve



Dental Chair Sleeve




Tubing Sleeve



Tray/Head Rest Sleeve



X-Ray Sensor Sleeve



Barrier Film Dispenser

CLINIC MISCELLANEOUS



Endo Stand



Endo Assist



Bar Stand 22-Bolt



Germicide Soaker



Setup Tray - Flat (without cover)



B-Lok Tray - Divided (without cover)



Implant Organizer



Instrument Mat



Mini Mat



Restorative Organizer



Composite Tab



Syringes Stand 20-Unit

Procedure Tab - Complete Set
(with all insert)

Procedure Tab



Tab Cap (with cover)



Sterile Endo Box

CLINIC MISCELLANEOUS



Compact Cassette



Compact Sterile Container



Rack & Tube



Lab Pass



Adjustable 6-Position Rack



Cabinet Tray



Wall Hanging Lab Pass



Tub Insert



Procedure Tub Cover



B-Lok Tray Cover



Tub Divider



Lab Pass Strip 18"

TEACHING AID



CONSUMABLE PRODUCTS

