

STERNMED®



Docking Station - Pumps

Docking Station





Docking Station - Pumps

Docking Station

SternMed offers docking stations for the horizontal pumps SPN 58 and IPN 56 in two different sizes. Up to 6 or 8 pumps can be integrated to one docking station. All pumps will have unified power.

The Sternmed docking stations are easy to install and designed to be convenient for operation and maintenance. Both versions are equipped with an infusion hook for support of different pumps at the same time.

POWERFUL FEATURES

- Easy to install, convenient for operation and maintenance
- Up to 8 pumps can be installed in one docking station
- Equipped with infusion hook, supports the different pumps at the same time
- Designed for integration of syringe and infusion pumps
- Balanced design with high stability and smooth movement
- One unified power cord for all pumps

TECHNICAL SPECIFICATIONS DOCKING STATION - PUMPS | Docking Station SPN 58, IPN 56

	6 Pumps	8 Pumps
Power	AC 110 - 230 V, 50 - 60 Hz	AC 110 - 230 V, 50 - 60 Hz
Height (mm)	1550 mm	1690 mm
Width (mm)	530 mm	530 mm
Depth (mm)	550 mm	550 mm
Weight (Kg)	20 kg	21 kg
Infusion hook capacity	4	4
Wheels	4 Anitstatic Double wheels with Brake	4 Anitstatic Double wheels with Brake
Unified Power	Yes	Yes