



ACEONE-ENDO



ACEONEDENT KOREA IND. CO.



Since our foundation in 1991,
We have grown up as a leading manufacturer and exporter in the field of Gutta Percha Points and Absorbent Paper Points in Korea.
We acquired ISO9001, ISO13485 and CE Mark Certification by SGS in United Kingdom and registered to FDA in USA.
“ACEONE-ENDO” brand, endodontic points is proved an excellent quality of products.

- 
- 2009** 04 Appointed as a "Geonggi Internet Trade Frontier Business"
- 2007** 04 Registered as a manufacturer of medical device to KGMP
12 Appointed as a "Geonggi Industry Family Business"
- 2005** 07 Appointed as a "Geonggi Internet Trade Frontier Business"
06 Appointed as a "Geonggi Industry Family Business"
05 ISO 13485 Certified
- 2003** 07 Selected as a "Promising exporter by Small Business Corporation"
- 1999** 11 ISO 9001/ISO 13488, CE MARK Certified
07 Registered trademark "ACEONE-ENDO" to Patent Bureau (No. 381895)
- 1996** 12 Listed to FDA, USA
- 1994** 05 Registered as a manufacturer of medical device to KFDA
- 1993** 03 Company name changed to ACEONEDENT KOREA INDUSTRIAL COMPANY
- 1991** 02 Established C&L Dental Supply Trading Co.


GUTTA PERCHA POINTS

Gutta Percha Points is being used as a filling material (permanent use) of the cleaned root canals of tooth intended to fill the root canal of a tooth manufactured in strictly accordance with internationally recognized standards.

- True to size and shape with excellent workability
- Stiff but flexible to provide ideal workability
- Long shelf life and freshness assured
- Hand-rolled from finest Gutta Percha
- Uniformed and non-distorting crimp




Coloring / Package Type



Color Coded and Length Marked



Square Vials



Endo-Aide (Round Type)

ISO Standard Size

ISO Size	Color	Taper	Coloring	Package Type	Points / Box
No. 8	Gray	Taper 0.02	<ul style="list-style-type: none"> · Color Coded · Color Coded and Length Marked 	<ul style="list-style-type: none"> · Slide Box (6Compartment) · Square Vials 	<ul style="list-style-type: none"> · Taper 0.02 Color Coded / Length Marked : 120 points · Taper 0.04 & 0.06 Color Coded : 60 points
No. 10	Purple				
No. 15	White				
No. 20	Yellow				
No. 25	Red	Taper 0.02, 0.04 & 0.06			
No. 30	Blue				
No. 35	Green				
No. 40	Black				
No. 45	White				
No. 50	Yellow				
No. 55	Red	Taper 0.02			
No. 60	Blue				
No. 70	Green				
No. 80	Black				
No. 90	White				
No. 100	Yellow				
No. 110	Red	Taper 0.02, 0.04 & 0.06			
No. 120	Blue				
No. 130	Green				
No. 140	Black				
Assorted No. 15~40	-				
Assorted No. 45~80	-		Taper 0.02		
Assorted No. 15~80	-	Taper 0.02	Color Coded	Endo-Aide (Round type)	240 points



Taper 0.02



Taper 0.04



Taper 0.06

Non - Standard Size -Accessory

ISO Size	Coloring	Package Type	Points / Box
X-Fine	None	Slide Box(6Compartment)	100 Points
Fine-Fine			
Medium-Fine			
Fine			
Fine-Medium			
Medium			
Medium-Large			
Large			
X-Large			
Assorted (XF, FF, MF, F, FM, M)			



ABSORBENT PAPER POINTS


Absorbent Paper Points is a device made from highly absorbent natural paper intended to clean the root canal of a tooth by insertion and removal supplied in sterilized condition.

- True to size and shape with excellent workability
- Uniformly hand-rolled to exact compliance to allowable tolerance
- Unsurpassed consistency in quality
- Rigid enough to be easily inserted into the moist root canal without curing or bending and be highly absorbent excellent radio-sterilized

Coloring / Package Type



Cell Packs




Color Coded and Length Marked




Endo-Aide (Round Type)

ISO Standard Size

ISO Size	Color	Taper	Coloring	Package Type	Points / Box
No. 8	Gray	Taper 0.02	<ul style="list-style-type: none"> · Color Coded · Color Coded and Length Marked 	<ul style="list-style-type: none"> · Slide Box (6Compartment) · Square Vials 	<ul style="list-style-type: none"> · Taper 0.02 Color Coded / Length Marked : 200 points · Taper 0.04 & 0.06 Color Coded : 60 points
No. 10	Purple				
No. 15	White				
No. 20	Yellow				
No. 25	Red	Taper 0.02, 0.04 & 0.06			
No. 30	Blue				
No. 35	Green				
No. 40	Black				
No. 45	White				
No. 50	Yellow				
No. 55	Red	Taper 0.02			
No. 60	Blue				
No. 70	Green				
No. 80	Black				
No. 90	White				
No. 100	Yellow				
No. 110	Red				
No. 120	Blue				
No. 130	Green	Taper 0.02, 0.04 & 0.06			
No. 140	Black				
Assorted No. 15~40	-	Taper 0.02			
Assorted No. 45~80	-	Taper 0.02			
Assorted No. 15~80	-	Taper 0.02	Color Coded	Endo-Aide (Round type)	400 points



Taper 0.02



Taper 0.04



Taper 0.06

Non - Standard Size -Accessory

ISO Size	Coloring	Package Type	Points / Box
X-Fine	None	6Compartment	200 Points
Fine			
Medium			
Coatse			
X-Coatse			

